


# CombiScreen® PLUS



- Urine test strips with ascorbic acid protected glucose and blood pad
- Perform highest quality standards
- Excellent price-performance ratio
- Automated and visual reading



combi <sup>®</sup> screen PLUS			Glucose	Ascorbic Acid	Ketones	Protein	pH-Value	Blood	Nitrite	Leucocytes	Spec. Gravity	Bilirubin	Urobilinogen
<b>Name</b>	<b>Art.-No.</b>	<b>Cont.</b>											
<b>Glucose PLUS</b>	94 501	50	■										
	94 101	100	■										
<b>Nitrit PLUS</b>	94 506	50							■				
<b>3 PLUS</b>	94 508	50	■			■	■						
	94 108	100	■			■	■						
<b>5+N PLUS</b>	94 535	50	■	■		■	■	■					
	94 135	100	■	■		■	■	■					
<b>5+Leuko PLUS</b>	94 517	50	■			■	■	■	■				■
	94 117	100	■			■	■	■	■				■
<b>9 PLUS</b>	94 115	100	■	■	■	■	■	■				■	■
	94 250	50	■	■	■	■	■	■				■	■
<b>9+Leuko PLUS</b>	94 200	100	■	■	■	■	■	■				■	■
	94 120	100	■	■	■	■	■	■				■	■
<b>10 SL PLUS</b>	94 120	100	■	■	■	■	■	■	■	■	■	■	■
<b>5 SYS<sup>2)</sup> PLUS</b>	94 109	100	■			■	■	■	■				
<b>7 SYS<sup>2)</sup> PLUS</b>	94 110	100	■			■	■	■	■				
<b>11 SYS<sup>1)</sup> PLUS</b>	94 100	100	■	■	■	■	■	■	■	■	■	■	■
<b>VET11 PLUS</b>	94 100V	100	■	■	■	■	■	■	■	■	■	■	■

<sup>1)</sup> System test strips for CombiScan<sup>®</sup> und Urilyzer<sup>®</sup> instruments  
<sup>2)</sup> System test strips for Urilyzer<sup>®</sup> instruments

**Ascorbic acid (Vitamin C)** is being ingested through food and nutritional supplements and passed unaltered in the urine. For this reason a high percentage

combi <sup>®</sup> screen classic line (without ascorbic acid prot.)			Glucose	Ascorbic Acid	Ketones	Protein	pH-Value	Blood	Nitrite	Leucocytes	Spec. Gravity	Bilirubin	Urobilinogen
<b>Name</b>	<b>Art.-No.</b>	<b>Cont.</b>											
<b>GP</b>	93 104	100	■			■							
<b>3</b>	93 108	100	■			■	■						
<b>GAK</b>	93 107	100	■	■	■								
<b>GPK</b>	93 105	100	■		■	■							
<b>5 B</b>	93 134	100	■		■	■	■	■					
<b>10 SL</b>	93 120	100	■		■	■	■	■	■	■	■	■	■
<b>11 SYS</b>	93 100	100	■	■	■	■	■	■	■	■	■	■	■

of the population shows a clearly increased ascorbic acid level in the urine. Depending on the percentage of ascorbic acid, lowered or even false negative results on the blood and glucose pad may occur. Due to that fact, necessary diagnostic or therapeutic actions might be delayed. In order to obtain precise results CombiScreen PLUS is equipped with an integrated ascorbic acid protection system.



#### Kidney and urinary tract diseases

Nitrite, Protein, Specific Gravity, Blood, Leucocytes and pH-Value.




#### Diabetes mellitus

Glucose and Ketones.




#### Liver and bile diseases

Bilirubin and Urobilinogen.



CombiScreen<sup>®</sup> Control PN: simple and safe control of the analyzers and test strips.



Urine analyzers  
 Urilyzer<sup>®</sup> 100, CombiScan<sup>®</sup> 500 and XL – the ideal choice for doctors' surgeries, hospitals and laboratories.

Distributor information

**Analyticon  
 Biotechnologies AG**

Am Muehlenberg 10  
 35104 Lichtenfels - Germany

Fon +49 (0) 6454 7991-0  
 Fax +49 (0) 6454 7991-30

info@analyticon-diagnostics.com  
 www.analyticon-diagnostics.com